

**CATALOGUE #:** 3IN5

**PRODUCT NAME:** Monoclonal anti-influenza virus type A (nucleoprotein)

**Recombinant  
MAbs:**

**FA32, FA35, FA38, FA58**

Recombinant antibody expressed in a mammalian cell line. Full-size IgG sequence derived from rabbit B cells.

**FA52**

Recombinant chimeric antibody expressed in a mammalian cell line. Composed of original wild type variable domains of rat derived MAb and human IgG1 constant domains.

**FA91, FA94**

Recombinant chimeric antibody expressed in a mammalian cell line. Composed of original wild type variable domains of sheep derived MAb and human IgG1 constant domains.

**MAbs *in vitro*:**

**FA17**

Mouse monoclonal antibody produced in bioreactor. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.

**MAbs *in vivo*:**

**F8, InA108, InA180, InA224, InA245**

Mouse monoclonal antibody produced in ascites. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.

**Immunogens:**

Recombinant nucleoprotein of influenza virus type A strain H1N1 A/California/07/2009 for FA17, FA32, FA35, FA38, FA52, FA58, FA91, FA94

Purified influenza virus type A strain H1N1 for F8, InA108, InA180, InA224, InA245

**Specificity:**

Influenza virus type A (nucleoprotein). No cross-reactivity to NP of influenza B virus.

The following strains were tested with FA17, FA32, FA35, FA38, FA52, FA58, FA91, and FA94 and they all reacted with these strains:

A/California/07/2009(H1N1)	A/Texas/50/2012(H3N2)
A/Taiwan/1/1986(H1N1)	A/Brisbane/10/2007(H3N2)
A/Beijing/262/1995(H1N1)	A/Singapore/1/1957(H2N2)
A/New Caledonia/20/1999(H1N1)	A/Tern/South Africa/1961 H5N3)
A/Solomon Islands/03/2006(H1N1)	A/Mexico/InDRE7218/2012(H7N3)
A/Hong Kong/45/2019(H3N2)	A/chicken/Nakorn-Patom/Thailand/CU-K2/2004(H5N1)
A/Panama/2007/1999(H3N2)	A/chicken/HongKong/NT142/2003(H9N2)
A/Wisconsin/67/2005(H3N2)	A/Anhui/1/2013(H7N9)

Testing was carried out using corresponding recombinant nucleoproteins or lysates of purified viral preparations.

All antibodies are not cross-reactive to influenza B virus (<0.1%). Testing was carried out using lysates of purified viral preparations of influenza B/Colorado/06/2017.

FA17, FA32, FA35, FA38, FA52, FA58, FA91, FA94 were tested with SARS-CoV-2 nucleoprotein and demonstrated no cross-reaction (<0.05%).

**CATALOGUE #:** 3IN5

**PRODUCT NAME:** Monoclonal anti-influenza virus type A (nucleoprotein)

**MAb isotypes:** IgG for FA32, FA35, FA38, FA58  
IgG1 for FA17, FA52, FA91, FA94, InA108, InA224  
IgG2a for F8  
IgG2b for InA245  
IgG3 for InA180

**Applications:** F8 can be used in immunocytochemistry. F8 inhibits viral reproduction after fatty acid acylation.  
InA108 and InA245 detect influenza A nucleoprotein in Western blotting.  
MAbs are working in ELISA.  
Recommended pairs for influenza A nucleoprotein detection:

Capture – Detection	Lateral flow (gold nanoparticles for detection)	Sandwich immunoassay
FA35 – FA17	+	+
FA52 – FA17	+	+
FA32 – FA17	+	+
FA38 – FA17	+	+
FA94 – FA17	+	+
FA58 – FA17	+	+
FA91 – FA17	+	+
InA108 – InA245		+
InA224 – InA245		+

**Purification:** Protein A chromatography for FA17, FA32, FA35, FA38, FA52, FA58, FA91, FA94, InA108, InA180, InA224, InA245  
Protein G chromatography for F8

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>) for FA17, FA32, FA35, FA38, FA52, FA58, FA91, FA94, F8, InA108, InA224, InA245  
50 mM sodium citrate, 150 mM NaCl, pH 6.0, 0.09 % azide (NaN<sub>3</sub>) for InA180

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.