

Datasheet

Blood coagulation and Anemia • Bone Metabolism • Cardiac Markers • Fertility and Pregnancy
Gangliosides • Hormone Markers • Immunology and Serology • Infectious Diseases • Inflammation
Kidney Diseases • Metabolic Syndrome • Microbial and Plant Toxins • Miscellaneous • Neuroscience
Thyroid Diseases • Tumor Markers • Veterinary

CATALOGUE #: 8CA72

PRODUCT NAME: CA72-4

Source: Human metastatic liver carcinoma.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies and HCV.

Applications: Immunogen, calibrator in immunoassays, in vitro studies.

Presentation: PBS, pH 7.4, 0.05 % of sodium azide (NaN₃).

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.
Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.