

CATALOGUE #: 8M79

PRODUCT NAME: Myelin, basic protein

Source:	Human brain. Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.
Analysis:	Purity > 95 % (SDS-PAGE)
Immuno reactivity:	Immunoreactivity is checked using poly- and monoclonal antibodies against MBP.
Purification:	Purified according to G. E. Deibler, Preparative Biochemistry 2(2) 139-165 (1972).
Presentation:	Lyophilized solid, essential salt free
Storage:	-20 °C (-15 ... -30 °C allowed)
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material.