

CATALOGUE #: 4I52

PRODUCT NAME: Monoclonal mouse anti-insulin-like growth factor binding protein-1 (IGFBP-1 (pp12))

MAbs <i>in vivo</i>:	G2 Mouse monoclonal antibody produced in ascites. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.
Immunogen:	IGFBP-1 purified from human amniotic fluid
MAb isotypes:	IgG2a for G2
Applications:	Indirect ELISA and Western blotting
Purification:	Protein G chromatography
Presentation:	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	<p>This product is sold for research or further manufacturing use only. Standard Laboratory Practices should be followed when handling this material.</p> <p>Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.</p>