

CATALOGUE #: 8HC51

PRODUCT NAME: HCV, NS4 antigen, recombinant

Source:	Expressed in <i>E. coli</i>
Description:	Recombinant antigen Contains C-terminal 10His-Cys-tag
Applications:	Immunoassays
Analysis:	Purity >90% Concentration is determined spectrophotometrically
Purification:	Metal affinity chromatography
Presentation:	50 mM Tris-HCl, pH 8.0, 2.5 mM EDTA, 5 % sucrose, 0.05 % SDS, 0.09 % sodium azide (NaN ₃)
Storage:	-70 °C (-65 ... -80 °C allowed)
Other information:	Avoid repeated freezing and thawing. It is recommended to aliquot the product after the first thawing.
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.