

CATALOGUE #: 8HC53

PRODUCT NAME: HCV, NS3 antigen, recombinant

Description: Recombinant antigen expressed in *E. coli*. Contains C-terminal 10His-Cys-tag.

Applications: Immunoassays

Analysis: Purity >92%

Concentration is determined spectrophotometrically at A280

Purification: Metal affinity chromatography

Presentation: 50 mM Tris, pH 8.0, 2.5 mM EDTA, 5 % sucrose, 0.09 % sodium azide (NaN₃)

Storage: -70 °C (-65 ... -80 °C allowed)

Other information: Avoid repeated freezing and thawing. It is recommended to aliquot the product after the first thawing.

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.