

CATALOGUE #: 502

PRODUCT NAME: Monoclonal mouse anti-ovine IgG

MAbs <i>in vivo</i>:	9E2 Mouse monoclonal antibody produced in ascites. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.
Immunogen:	Purified ovine IgG
Specificity:	Ovine IgG Cross-reacts with bovine, capricorn and deer IgG.
MAb isotypes:	IgG1 for 9E2
Applications:	Detection of specific “late” antibodies by indirect ELISA.
Purification:	Protein A chromatography
Presentation:	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	This product is sold for research or further manufacturing use only. Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.