

**CATALOGUE #:** 2FSH2

**PRODUCT NAME:** Monoclonal mouse anti-Follicle stimulating hormone (FSH), beta chain

**MAbs *in vivo*:** **F2**  
Mouse monoclonal antibody produced in ascites. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.

**Immunogen:** Purified natural FSH beta

**Specificity:** High affinity for the beta chain of FSH.  
No cross-reactivity with LH, TSH or HCG.

**MAb isotypes:** **IgG1** for F2

**Applications:** EIA, Western blotting  
Recommended pair for sandwich immunoassay:

Capture	Detection
F2	F1cc

(F1cc is under Cat.# 2H8)

**Purification:** Ion exchange chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.