

CATALOGUE #: 6L56

PRODUCT NAME: Monoclonal mouse anti-CD56

MAb:	LT56cc Mouse monoclonal antibody produced in bioreactor. Hybridoma clone derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice.
Immunogen:	KG1a cell line
Specificity:	Surface cell antigen CD56 (MW 175-220 kDa).
MAb isotype:	IgG2a for LT56cc
Presentation:	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	<p>This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.</p> <p>Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.</p>