

CATALOGUE #: 8G16-5b

PRODUCT NAME: Monosialoganglioside GM4, bovine

Source: Bovine brain tissue.

Application: Ganglioside converted to the galactosylceramide by bacterial or other sialidase.

Analysis: Purity more than 98 % by HPTLC. MW 1017

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.