



# DataSheet

**CATALOGUE #:** 2GH1

**PRODUCT NAME:** Monoclonal mouse anti-ghrelin

---

**MAbs:** Ghr5, Ghr6

**MAb isotypes:** **IgG2a** for MAb Ghr5  
**IgG2b** for MAb Ghr6

**Applications:** MAbs are working in ELISA.

**Purification:** Protein G chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.