



DataSheet

CATALOGUE #: 2LTB2

PRODUCT NAME: Monoclonal mouse anti- *Escherichia coli* heat-labile enterotoxin B-chain

MAbs: BB2, BG12

Specificity: *Escherichia coli* heat-labile enterotoxin B-chain

MAb isotypes: **IgG1** for MAb BG12

IgG2b for MAb BB2

Applications: EIA and Western blotting

Purification: Salt precipitation. Purity \geq 80 % (SDS-PAGE).

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN_3)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi