



DataSheet

CATALOGUE #: 2S6

PRODUCT NAME: Monoclonal mouse anti-*Staphylococcus aureus* enterotoxin type G

MAbs: SEG-59, SEG-16

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with *S. aureus* enterotoxin type G.

Specificity: *S. aureus* enterotoxin type G

MAB isotypes: **IgG1** for MAb SEG-59
IgG2a for MAb SEG-16

Applications: Recommended pair for sandwich immunoassay (capture-detection):
SEG-59 – SEG-16
Both MAbs are working in Western blotting.

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi