



# DataSheet

**CATALOGUE #:** 2S7

**PRODUCT NAME:** Monoclonal mouse anti-*Staphylococcus aureus* enterotoxin type A

---

**MAbs:** C4, E11, H5, H10

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with *S. aureus* enterotoxin type A.

**Specificity:** *S. aureus* enterotoxin type A

**MAB isotypes:** **IgG1** for MAbs C4, H5

**IgG2a** for MAb E11

**IgG2b** for MAb H10

**Applications:** ELISA and Western blotting.

Recommended pairs for sandwich immunoassay (capture-detection):

H5 - C4

H10 - C4

C4 - H5

E11 - C4

**Purification:** Protein G chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)