



DataSheet

CATALOGUE #: 3V1

PRODUCT NAME: Monoclonal mouse anti-vaccinia virus

MAbs: TV43, TV46

Specificity: MAbs are cross-reacting with smallpox virus.

MAB isotypes: **IgG2a** for MAbs TV43 and TV46

Applications: MAbs are working in ELISA and Western blotting.

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi