



DataSheet

CATALOGUE #: 4CA19

PRODUCT NAME: Monoclonal mouse anti-CA19-9

MAbs: X91

Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with affinity purified human CA19-9 antigen.

Specificity: CA19-9 antigen

MAB isotypes: **IgM** for MAb X91

Applications: Detection of CA19-9 antigen in immunoassays. MAb is working in immunohistochemistry.

Purification: Ion exchange chromatography

Presentation: 0.15 M NaCl, 10 mM Tris, pH 7.5, 0.05 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi