



DataSheet

CATALOGUE #: 8G16-10b

PRODUCT NAME: Trisialoganglioside GT1b, bovine

Source: Bovine brain tissue.

Applications: Useful as an inhibitor of Con A-stimulated mutagenesis in murine T-cells. Marker of lymphocytes from cancer patient. Readily convert to the GM1 ganglioside by using bacterial or other sialidase from mammalian tissue without detergent.

Analysis: Purity more than 98 % by HPTLC. MW 2109

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi