



DataSheet

CATALOGUE #: 8G16-19b

PRODUCT NAME: Ganglioside GM1 sphinganine, bovine

Source: Bovine brain tissue.

Analysis: Purity more than 98 % by HPTLC. MW 1548
Contains C18:0, C20:0 sphingosine only and saturated fatty acid residues (mostly palmitoyl)

Presentation: Lyophilized from chloroform-methanol solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi