



DataSheet

CATALOGUE #: 4IG8

PRODUCT NAME: Monoclonal mouse anti-insulin-like growth factor binding protein-1 (IGFBP-1)

MAbs: G5F8, C7B9

Specificity: MAbs do not cross-react with human insulin or proinsulin

MAb isotypes: IgG1 for MAbs C7B9, G5F8

Applications: Antibodies are working in ELISA and Western blotting.
Recommended pair for sandwich ELISA (capture-detection):
C7B9 – G5F8

Purification: Protein A chromatography

Presentation: PBS, pH 7.4, 0.04 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.