



DataSheet

CATALOGUE #: 3C13

PRODUCT NAME: Monoclonal mouse anti-cyclosporin A

MAbs: CSZ22

Specificity: Cyclosporin A
MAb CSZ22 recognizes cyclosporin A in buffer solutions and whole blood extracts.

MAb isotypes: **IgG1** for MAb CSZ22

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi