



# DataSheet

**CATALOGUE #:** 3L1

**PRODUCT NAME:** Monoclonal mouse anti-*Listeria monocytogenes*

---

**MAbs:** LZF7, LZH1

**Specificity:** MAbs recognize *Listeria monocytogenes* outer membrane fraction and intact *L. monocytogenes* cells.

**MAb isotypes:** **IgG1** for MAb LZH1  
**IgG2a** for MAb LZF7

**Applications:** Indirect ELISA and Western blotting  
Recommended pair for sandwich immunoassay (capture-detection): LZH1 – LZF7.

**Purification:** Protein G chromatography.

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)