



DataSheet

CATALOGUE #: 4GP31

PRODUCT NAME: Monoclonal mouse anti-glycogen phosphorylase isoenzyme BB (GPBB)

MAbs: 1G6

Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with human GPBB.

Specificity: BB isoenzyme

MAb isotypes: IgG2b for MAb 1G6

Applications: GPBB immunoassay, immunoaffinity purification, Western blotting.

Purification: Protein A chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi