



DataSheet

CATALOGUE #: RC4TC2

PRODUCT NAME: Recombinant anti-cardiac troponin complex

MAbs: RecChim20C6

Recombinant full-sized chimeric 20C6 expressed in mammalian cells. Original wild type 20C6 variable domains and human IgG1 constant domains are combined in this construct.

Specificity: Human cardiac troponin complex

Applications: Human cardiac troponin complex sandwich immunoassays

Purification: Protein A chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi