



DataSheet

CATALOGUE #: 8C10b

PRODUCT NAME: Calmodulin, bovine

Source: Bovine brain tissue

Application: Biochemical and immunochemical investigations.

Analysis: Purity > 95% (native PAGE in reducing conditions with EGTA).

Purification: Ion-exchange chromatography followed by gel-filtration.

Presentation: Lyophilized from protein solution with 2 mM EDTA.

It is recommended to reconstitute this product with deionized water to its initial concentration.

Storage: Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi