



# DataSheet

**CATALOGUE #:** 8C10h

**PRODUCT NAME:** Calmodulin, human

---

**Source:** Human brain tissue. Blood samples from tissue donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

**Application:** Biochemical and immunochemical investigations.

**Analysis:** Purity > 95 % (native PAGE in reducing conditions with EGTA).

**Purification:** Ion-exchange chromatography followed by gel-filtration.

**Presentation:** Lyophilized from protein solution with 2 mM EDTA.  
It is recommended to reconstitute this product with deionized water to its initial concentration.

**Storage:** Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)