



DataSheet

CATALOGUE #: 8CA15

PRODUCT NAME: CA15-3, standard grade

Source: Human milk.
Blood sample from the donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Applications: Calibrator in immunoassays, in vitro studies.

Purification: High-speed centrifugation, salt precipitation and delipidization.

Presentation: Solution in 1 M NaCl, containing 0.09 % of sodium azide as preservative.

Storage: -20°C (-15 ... -30 °C allowed)

Other information: Avoid repeated freezing and thawing. It is recommended to aliquot the product after the first thawing.

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.
Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligent 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi