



DataSheet

CATALOGUE #: 8CA19

PRODUCT NAME: CA19-9

Source: Human metastatic liver carcinoma.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Applications: Immunogen, calibrator in immunoassays, in vitro studies.

Specific activity: Specific activity confirmed by enzyme immunoassay.

Presentation: PBS, pH 7.4, 3 % sucrose, 0.05 % of sodium azide (NaN₃).

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.
Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi