



DataSheet

CATALOGUE #: 8CT8

PRODUCT NAME: Canine thyroglobulin

Source: Canine thyroid gland

Applications: To be used as an antigen in immunoassays for determination of thyroglobulin autoantibodies in canine serum.

Analysis: Purity > 90% (SDS-PAGE)
Immunological identity confirmed by reaction with monoclonal antibody that is specific for the canine thyroglobulin.

Presentation: Lyophilized from PBS, pH 7.4
It is recommended to reconstitute this product with deionized water to its initial concentration.

Storage: Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)

Other information: Avoid repeated freezing and thawing. It is recommended to aliquot the product after reconstitution.
It is recommended to use PBS, pH 7.4, to prepare diluted samples.

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi