



# DataSheet

**CATALOGUE #:** 8CTS5

**PRODUCT NAME:** Recombinant canine thyroid-stimulating hormone (TSH)

---

**Source:** Expressed in a mammalian cell line

**Description:** Canine TSH is a heterodimeric glycoprotein consisting of two non-covalently linked subunits, an alpha subunit (96 amino acid residues, aar) and a beta subunit (118 aar).

Recombinant canine TSH is produced by co-expression of the alpha (UniProt Q9XSW8) and beta (UniProt P54828) subunits of TSH. Beta subunit contains alanine instead of valine at position 81 (Val81Ala polymorphism in canine TSH beta subunit has been reported).

**Applications:** Calibrator or standard in immunoassays. Biological activity has not been determined.

**Analysis:** Purity > 90% (SDS-PAGE in reducing conditions)

Product may contain free beta subunit

Concentration is determined by Lowry method using BSA as a standard

**Purification:** Affinity chromatography

**Presentation:** Lyophilized from 10 mM K-phosphate, pH 7.4, 150 mM NaCl, 100 mM D-mannitol, 0.1% CHAPS.

It is recommended to reconstitute this product with deionized water to its initial concentration.

**Storage:** Lyophilized -20°C (-15 ... -30 °C allowed)

Reconstituted -70°C (-65 ... -80 °C allowed)

**Other information:** Avoid repeated freezing and thawing.

It is recommended to use 10 mM K-phosphate, pH 7.4, 150 mM NaCl, 100 mM D-mannitol, 0.1% CHAPS, to prepare diluted samples.

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligence 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)