



DataSheet

CATALOGUE #: 8M79

PRODUCT NAME: Myelin, basic protein

Source: Human brain.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Analysis: Purity > 95 % (SDS-PAGE)

Immuno reactivity: Immunoreactivity is checked using poly- and monoclonal antibodies against MBP.

Purification: Purified according to G. E. Deibler, Preparative Biochemistry 2(2) 139-165 (1972).

Presentation: Lyophilized solid, essential salt free

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.