



DataSheet

CATALOGUE #: 8RP7

PRODUCT NAME: Retinol-binding protein 4 (RBP4) from human plasma, complexed with prealbumin

Source:	Pooled human plasma Donor blood samples were tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV
Composition:	Human RBP4 in complex with human prealbumin (transthyretin)
Applications:	Calibrator in immunoassays
Analysis:	Purity > 70 % (SDS-PAGE)
Immunological activity:	Immunoreactivity is checked using monoclonal antibodies specific to RBP4
Purification:	Affinity chromatography and size-exclusion chromatography
Presentation:	Lyophilized from phosphate-buffered saline (10 mM K-phosphate, pH 7.4, 150 mM NaCl) It is recommended to reconstitute this product with deionized water to its initial concentration.
Storage:	Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi