



# DataSheet

**CATALOGUE #:** 8G16-1b

**PRODUCT NAME:** Asialoganglioside GM1, bovine

---

**Source:** Bovine brain tissue.

**Analysis:** Purity more than 98 % by HPTLC. MW 1263

**Presentation:** Lyophilized from chloroform-methanol solution.

**Storage:** Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.