



DataSheet

CATALOGUE #: 8G16-12b

PRODUCT NAME: Tetrasialoganglioside GQ1b, bovine

Source: Bovine brain tissue.

Applications: A typical marker of the peripheral area of biopsy specimens for human glyomas. Can be used as high-sialylated substrate for neuraminidase.

Analysis: Purity more than 98 % by HPTLC. MW 2391

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi