



DataSheet

CATALOGUE #: 8G16-17b

PRODUCT NAME: Disialoganglioside GD1a-NAcGal, bovine

Source: Bovine brain tissue.

Applications: One of the minor neuronal gangliosides. Has the terminal GalNAc-residue with beta linkage to the terminal Gal-residue of GD1a. It is stable against the majority of sialidases, including the sialidase from *Vibrio cholerae*.

Analysis: Purity more than 98 % by HPTLC. MW 2030

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi