



# DataSheet

**CATALOGUE #:** 3BA18

**PRODUCT NAME:** IgG fraction of rabbit serum to spore antigen of *Bacillus anthracis*

---

**Immunogen:** Spore extract (MW 90 kD)

**Host animal:** Rabbit

**Specificity:** Spore antigen of *Bacillus anthracis*. Cross-reactivity with extracts from spore and vegetative cells of other *Bacillus* species is less than 5 %.

**Applications:** Detection of *Bacillus anthracis* spores

**Purification:** Protein A chromatography

**Presentation:** PBS, pH 7.2, 0.005 % merthiolate and 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)