



# DataSheet

**CATALOGUE #:** PPC3

**PRODUCT NAME:** Polyclonal anti-procalcitonin (PCT)

---

**Immunogen:** Human recombinant full-length procalcitonin

**Host animal:** Goat

**Specificity:** Specific to human PCT. Tested with human recombinant full length PCT. Cross-reactivity with other possible proteins has not been tested.

**Applications:** Human procalcitonin immunoassay

**Purification:** Immunoaffinity method utilizing PCT-sepharose as an affinity matrix.

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligence 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)