



DataSheet

CATALOGUE #: PRP4

PRODUCT NAME: Polyclonal anti-canine C-reactive protein (cCRP)

Description: IgG fraction of goat antibodies

Immunogen: Canine recombinant CRP

Host animal: Goat

Specificity: Specific to canine CRP, no cross-reactivity to other serum proteins in sandwich immunoassay.

Applications: Canine CRP immunoassays

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi