



# DataSheet

**CATALOGUE #:** 8BFP

**PRODUCT NAME:** BNP and NT-proBNP free plasma, EDTA

---

**Source:** Pooled normal human EDTA plasma.  
Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies and HCV.

**Application:** BNP and NT-proBNP free plasma is free of whole proBNP molecule and any proBNP fragments (BNP and NT-proBNP). It can be used for preparation of any proBNP-related standards and calibrators for immunoassays based on detection of any kinds of proBNP-related peptides.

**Purification:** Immunoaffinity chromatography of pooled normal human plasma.  
Gel for affinity column was prepared by coupling anti-BNP and anti-NT-proBNP monoclonal antibodies with different epitope specificity with BrCN-activated Sepharose CL 4B.

**Presentation:** Frozen liquid.

**Storage:** -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)