



# DataSheet

**CATALOGUE #:** 8FFS

**PRODUCT NAME:** Fatty acid binding protein (FABP) free serum

---

**Source:** Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

**Application:** FABP free serum can be used for preparation of human FABP standards and calibrators.

**Analysis:** According to ELISA testing FABP free serum does not contain more than 0.5 ng/ml of human FABP.

**Purification:** Immunoaffinity chromatography.  
Gel for affinity column was prepared by coupling anti-FABP monoclonal antibodies (several monoclonal antibodies specific for different epitopes of FABP molecule) with BrCN-activated Sepharose CL 4B.

**Presentation:** Frozen liquid.

**Storage:** -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligent 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)