



DataSheet

CATALOGUE #: 8MFS

PRODUCT NAME: Myoglobin free serum

Source: Pooled normal human serum.
Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Application: Myoglobin free serum can be used for preparation of human myoglobin standards and calibrators.

Analysis: According to ELISA testing myoglobin free serum does not contain more than 1 ng/ml of human myoglobin.

Purification: Immunoaffinity chromatography of pooled normal human serum.
Gel for affinity column was prepared by coupling anti-myoglobin monoclonal antibodies (several monoclonal antibodies specific to the N-terminal, C-terminal and central parts of myoglobin molecule) with BrCN-activated Sepharose CL 4B.

Presentation: Frozen liquid

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi