



DataSheet

CATALOGUE #: 8MPFS

PRODUCT NAME: Myeloperoxidase free serum

Source: Pooled normal human serum.
Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Application: Myeloperoxidase free serum can be used for preparation of human myeloperoxidase standards and calibrators.

Analysis: No detectable amounts of human myeloperoxidase.

Purification: Affinity chromatography of pooled normal human serum.
Gel for affinity column was prepared by coupling anti-myeloperoxidase monoclonal antibodies (several monoclonal antibodies with different epitope specificity) with BrCN-activated Sepharose CL 4B.

Presentation: Frozen liquid

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi