



# DataSheet

**CATALOGUE #:** 3L15

**PRODUCT NAME:** Monoclonal mouse anti-*Legionella pneumophila* LPS

---

**MAbs:** 2F10, 5F4RC

Hybridoma clones have been derived from hybridization of mouse Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with membrane fraction of *Legionella pneumophila*, strain Philadelphia 1.

**Specificity:** MAbs react with LPS membrane of *Legionella pneumophila*. MAbs do not cross-react with *Legionella pneumophila* serogroups 2, 6, 7, 10 and 11, *Legionella bozemanii*, *Legionella lonbeachae*, *Legionella micdadei*, *Corynebacterium diphtheria*, *Pseudomonas fluorescens*, *Pseudomonas cepacia*, *Pseudomonas aeruginosa*, *Bordetella pertusis*, *Leptospira interrogans* (Australia, Pomona, Icterogemorrhagia), *Toxoplasma gondii*, *Hemophilus influenza (type B)*, *Hemophilus ducreyi*, *Brucella abortus*, *Bacillus anthracis*, *Klebsiella pneumoniae*, *Listeria monocytogenes*, *Mycobacterium tuberculosis*, *Neisseria meningitides*, *Pasteurella multocida*, *Proteus vulgaris*, *Shigella sonnei*, *Staphylococcus faecalis*, *Streptococcus puogenes*, *Streptococcus pneumoniae*, *Yersinia pseudotuberculosis*, *Yersinia pestis EV76*, *Salmonella typhi*, *Escherichia coli* K88 and *Francisella tularensis* 15.

Cross-reaction with *Legionella pneumophila* serogroups 3 and 4 has not been tested.

**MAb isotypes:** IgG3 for MAbs 5F4RC and 2F10

**Applications:** Detection of *Legionella pneumophila* in ELISA.

**Purification:** Protein G chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)