



DataSheet

CATALOGUE #: 3AV13

PRODUCT NAME: Monoclonal mouse anti-adenovirus hexon

MAbs: 8C4

Specificity: Hexon antigen of human adenovirus.

MAb isotypes: IgG2a for MAb 8C4

Applications: Detection of human adenovirus. MAb 8C4 can be used in ELISA, immunodiffusion, immunohistochemistry and lateral flow.

Purification: Protein A chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligent 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi